

# Gonad / ovarial-shields

## Flexible Gonad-Shields

For male patients

Optimal protection of the genitals. Allowing self application by the patient. Made from multiply lead rubber sheets with vinyl cover for easy cleaning. Lead equivalent: 1.0 mm Pb

Application during radioscropy of pelvis and especially of hip-joint, thigh and lumbar vertebral column.



871003 - M Set = 3 Sizes

sizes	colour
small below 9 years	scarletred
medium from 9-14 years	petrolgreen
large over 15 years	navyblue

871013-M Set = 3 sizes with hipster



871003-M/H

on request: hygiene pouches set for all sizes

## Flexible Ovarial-Shields

for female patients

Application during radioscropy of pelvis and especially of hip-joint, thigh and lumbar vertebral column.

Lead equivalent: 1,0 mm

871003 - W

Set = 3 sizes

sizes	colour
small below 9 years	scarletred
medium from 9-14 years	petrolgreen
large over 15 years	navyblue




## Ovarial-Protector OP

The special form guarantees protection of higher positioned Ovarians.

Pure lead covered with hygienic PVC and easy to clean.

Special size for babies and 5 standard sizes available - with engraved size and age group.



S Special size for infants colour: jeans

I below 4 years	colour: light blue
II from 4-8 years	colour: gray
III from 8-12 years	colour: green
IV from 12-16 years	colour: beige
V over 16 years	colour: silver gray

87005-W

Set = 5 sizes I - V